



TECHNICAL CONSIDERATIONS FOR AN OPTIMIZED JOINT PROTOCOL


Raúl San José Estépar, PhD
Brigham and Women's Hospital



Introduction: The compromise of a joint protocol



The nature of the challenge



Physics of the scanner

Scanner Configuration

Tube current (dose)

Table centering



Software of the scanner

Current Modulation

Noise

Reconstruction software

- Recon Kernel
- Iterative Approaches

**PSF +
Noise**

- Slice Thickness


PSF

Minimize Risk
Maximize Quality


Recon vs. Tube Current: What is critical?

- *Metrology*

Point Spread Function



Noise



Recon vs. Tube Current: What is critical?


- *Densitometry*

Point Spread Function


Noise

Reconstruction Technique


FBP: STANDARD Kernel




FBP: LUNG Kernel




Iterative: ASIR 70




Iterative: VEO







Motion: Two frequencies and one finite beam




Reasons for Failure



Cardiac Motion: Is Gating needed?



Agatston Protocol



Volumetric Protocol

Where are we now?: NLST Protocol

Scanner Configuration				
Scanner Make	GE		Siemens	
Scanner Model	Light Speed 16	VCT(64)	Sensation 16	Senstation 64
Rotation Time (s)	0.5	0.5		
Det. Configuration	16x1.25	64x0.625	16x0.75	32x0.6
Pitch	1.375	0.984	1.5	1
Speed (mm/rot)	27.5	39.37	18	19.2
kVp	120	120	120	120
mAs	40	25	75	50
Dose Modulation				
Reconstruction				
Algorithm	Standard	Standard	B30	B30
Thickness (mm)	2.5	2.5	2	2
Interval (mm)	2	2	1.8	1.8
DFOV (cm)				

Where are we now?: COPDGene Protocol

Scanner Configuration					
Scanner Make	GE		Siemens		
Scanner Model	Light Speed 16	Discovery HD750	Definition	Definition AS+	Definition Flash
Rotation Time (s)	0.5	0.5	0.5	0.5	0.5
Det. Configuration	16x0.625	64x0.625	64x0.6	128x0.6	128x0.6
Pitch	1.375	1.375	1.1	1	1
Speed (mm/rot)	13.75	13.75	13.2	38.4	38.4
kVp	120	120	120	120	120
mAs	200	200	200	200	200
Dose Modulation	Auto-mA Off	Off	CARE Dose 4D Off	CARE Dose 4D Off	CARE Dose 4D Off
Reconstruction					
Algorithm	Standard	Standard	B31f	B31f	B31f
Thickness (mm)	0.625	0.625	0.75	0.75	0.75
Interval (mm)	0.625	0.625	0.5	0.5	0.5
DFOV (cm)	Lungs	Lungs	Lungs	Lungs	Lungs

Where are we now?: CAC Protocol

Parameter	Scanner						
	Imatron (Electron-Beam)	LightSpeed Plus	MX8000	Volume Zoom	Volume Zoom	Aquilion	Sensation 64
Acquisition mode	Sequential	Sequential	Sequential	Sequential	Spiral	Sequential	Spiral
Electrocardiographic synchronization	Prospective	Prospective	Prospective	Prospective	Retrospective	Prospective	Retrospective
Peak voltage (kVp)	130	120	120	120	120	120	120
Rotation time (sec)	0.1 (scanning time)	0.5	0.5	0.5	0.5	0.5	0.33
Tube current–time product (mAs)*	63 (fixed)	25, 70, 145	10, 30, 65	20, 55, 135	20, 50, 115	20, 45, 90	20, 70, 145
Detector configuration (mm)	1 × 3	4 × 2.5	4 × 2.5	4 × 2.5	4 × 2.5	4 × 3	64 × 0.6
Section thickness (mm)	3	2.5	2.5	2.5	3	3	3
Table feed (millimeters per rotation)	3	10	10	10	3.75	12	3.84
Pitch	1	1	1	1	0.375	1	0.2
Reconstruction algorithm	Sharp	Standard	B	B35f	B35f	FC01	B35f
Effective dose (mSv) [†]	1.0	1.7	0.6	1.1	2.5	1.4	5.25

McCollough, Coronary Artery Calcium: A Multi- institutional, Multimanufacturer International Standard for Quantification at Cardiac CT Radiology, 2007

Emphysema: One marker for three protocols

LAA%-950


LAA%-925

LAA%-910

COPDGene

N=4059

Mean FEV1%:
97.40 ±11.51%




NLST

Dransfield et al,
Chest, 2007


N=359

Mean FEV1%:
74±18%



FHS Cardiac

N=3409



Median Value


	LAA%-950	LAA%-925	LAA%-910
COPDGene	0.95%	5.65%	13.32%
NLST	3.5%	12.6%	24.8%
FHS Cardiac	2.0%	7.1%	14.27%

The effect of multiple thresholds


COPDGene

N=4059

LAA%-950 vs LAA%-925





LAA%-950 vs LAA%-910





NLST

N=359



FHS Cardiac

N=3409




FHS: One attempt to a join protocol


	Heart Scan	Full Lung Scan
Manufacturer	GE Light Speed Ultra	GE Discovery STE
Voltage (kV)	120	120
Dosage (mAs)	104.32	68.1
Slice Width	2.5	0.625
Reconstruction	STANDARD	LUNG

Emphysema in Full Lung Scan


LAA%-950



LAA%-925





LAA%-910




N=2795

Median: 29.8% 35.91% 39.80%



Heart Scan vs Full Lung Scan with Gating



Conclusions

- Effects of reconstruction overpowers the effects of dose
 - Standardized reconstruction protocols that balances the trade-offs
 - Iterative reconstructions need to be validated
- Gating is an imperfect solution
- Quantitative markers need to be adapted to exploit invariances across scanning protocols and scanner generations