

Progress in Understanding Sources in Variance with Quantitative Imaging

Nicholas Petrick, Marios A. Gavrielides,
Rongping Zeng

Center for Devices and Radiological Health, U.S. Food and Drug
Administration

Purpose/Outline

- Discuss quantitative CT imaging research progress within OSEL/DIAM & through collaborative efforts
- Outline
 - Part 1: Overlapping reconstruction in CT volume estimation
 - Part 2: Update on “final” QIBA 1A reader study analysis

Part 1

Effect of overlapping CT reconstructions

RSNA 2011


Marios Gavrielides, OSEL/DIAM

Purpose

- To compare overlapping reconstruction to contiguous slice reconstruction
 - For the task of estimating nodule volume

Reconstructed CT slices

- Overlapping slice reconstruction



- Contiguous reconstruction



Reconstructed CT slices

- Overlapping reconstruction
 - Expressed as % of slice thickness
 - 50% overlap
- Contiguous reconstruction
 - Reconstruction interval = slice thickness
 - 0% overlap
- Same dose with each recon techniques

Trade-off between overlapping & contiguous recons

- Practical issues

- Overlap increases no. of slices
 - Reading time, storage size increases

- Improved z-dim resolution

- Overlap minimizes partial-volume effects
- Maximize contrast for central lesion slice*

*Kalender, JCAT, 1994., Wang, Med Phys, 1994, Leung J Thor Imag., 1997

Reconstructed CT slices


- Studies show improved detection with overlapping slices*
 - Effect on volume estimation not examined
- Compare overlapping to contiguous reconstructed CT slices for volume estimation task
 - CT data from our anthropomorphic phantom containing synthetic nodules

*Urban, AJR, 1993., Diederich, Eur Radiol 1999., Buckley, Radiol, 1995.

Imaged anthropomorphic phantom

- Limitations

- No airway, no lung parenchyma



Inserted synthetic nodules

- Spherical nodules


- Size: 5, 10, 20 mm
- Density: -630HU, +100HU
- Attached to lung vasculature



- CT Scanner

- Siemens 64-slice, Somatom 64
 - Images collected at the Mallinckrodt Institute of Radiology (MIR), Washington University School of Medicine

Collection protocol



Nodule Volume Estimator

- Matched filter-based approach*
 - Incorporates properties of imaging system

*Gavrielides et al., IEEE TMI, vol. 29, n. 10, pp. 1795-1807, 2010.

Analysis plan

- Comparison of biases

- Normalized size

- $Size_{Norm} = \frac{Size_{Est} - Size_{True}}{Size_{True}} \bullet 100\%$

- Relative bias

- $Bias_{Rel} = E\{Size_{Norm}\}$

Results

Function of size

Nodule Size mm	Overlapping Bias _{Rel} (Std)	Contiguous Bias _{Rel} (Std)
All sizes (N [*] =1080)	3.2 (15.0)	8.8 (29.2)
5 (N=360)	4.4 (24.1)	29.7 (33.2)
10 (N=360)	1.2 (7.7)	5.2 (9.9)
20 (N=360)	4.0 (5.6)	-8.4 (24.9)

- N= 1080
 - 3 sizes x 2 pitches x 3 exposures, 2 densities x 3 slice thicknesses x 10 repeats

Results

Fixed nodule to slice thickness ratio

Nodule Size/ Slice Thick Ratio	Overlapping Bias _{Rel} (Std)	Contiguous Bias _{Rel} (Std)
3.3 (N=240)	-1.2 (13.9)	7.3 (7.3)
6.6 (N=360)	0.5 (10.8)	-11.5 (23.3)
13.3 (N=240)	2.0 (7.5)	4.2 (8.5)

- 3.3: 5 mm @ 1.5 mm ST, 10mm @ 3 mm ST
- 6.7: 5 mm @ 0.75 mm ST, 10 mm @ 1.5 mm ST, 20mm @ 3 mm ST
- 13.3: 10 mm @ 0.75 mm ST, 20mm @ 1.5 mm ST

Part 1 summary

- Reduced Bias_{Rel} for overlapping slices compared to contiguous
 - Finding consistent across nodule sizes, nodule densities & imaging protocols
- Studies needed to determine whether increase in accuracy/precision outweigh workload issues
 - Contiguous for interpretations?
 - Overlap for nodule evaluation only?

Part 2

Update on QIBA 1A reader study analysis

QIBA volumetric CT Working Group

Purpose

- General purpose

- To provide supporting data and a methodological pathway for validating the technical performance of lesion sizing techniques

- Specific purpose

- To estimate the bias and variance of radiologists measuring the size of simple and complex synthetic nodules through a controlled reader study

Study design


- Data from QIBA 1A reader study
 - Readers size synthetic nodules from CT scan data
 - 10 synthetic nodules
 - Types
 - 10 mm Sphere
 - 20 mm Sphere
 - 20 mm Ellipsoid
 - 10 mm Lobulated
 - 10 mm Spiculated
 - Density
 - -10 HU
 - +100HU
- Readers sized nodules from 40 CT datasets
 - 10 nodules X 2 slice thickness X 2 repeat scans



Spherical



Ellipsoid




Lobulated



Spiculated

Study design



- Anthropomorphic thorax phantom
 - Kyotokagaku Incorporated, Tokyo, Japan

Reading protocol

- Readers

- 6 radiologists familiar with evaluating lesion response in drug trials

- Sizing Methods

- 1D technique (linear distance)
 - Largest in-slice diameter for the lesion
 - Based on RECIST criteria
- 2D technique (area)
 - Largest in-slice diameter for the lesion
 - Largest perpendicular diameter within same slice
 - Based on WHO criteria
- 3D technique (volume)
 - Semi-automated volumetric measurement tool

Analysis plan

- Comparison of biases and variances

- Normalized size

- $Size_{Norm} = \frac{Size_{Est} - Size_{True}}{Size_{True}} \bullet 100\%$


- Relative bias

- $Bias_{Rel} = E\{Size_{Norm}\}$


- Relative standard deviation

- $Std_{Rel} = Std\{Size_{Norm}\}$

Reading process




1D Size_{Norm} = -9.1%




2D Size_{Norm} = -44.6%

Reading process



3D In-slice



3D Size_{Norm} = +36.1%


Analysis plan

- ANOVA & goodness-of-fit statistics (R^2) to identify important contributing factors
- Multiple comparisons of difference in Bias_{Rel} among sizing methods
 - Bias_{Rel} & bias comparisons based on t-test within each subgroup
- Multiple comparisons of difference in relative variability among sizing methods
 - Std_{Rel} & Std comparisons based on bootstrapping within each subgroup

Results

ANOVA single factor


Source	Prob>F	Significant (p≤0.05)
Nodule Type	0	✓
Sizing Method	0	✓
Readers	0	✓
Nodule Density	0.0074	✓
Slice Thickness	0.0084	✓
Nodule Set	0.35	✗
Reading Session	0.7032	✗



Results

ANOVA single factor


Source	Prob>F	Significant (p≤0.05)
Nodule Type	0	✓
Sizing Method	0	✓
Readers	0	✓
Nodule Density	0.0074	✗
Slice Thickness	0.0084	✗
Nodule Set	0.35	✗
Reading Session	0.7032	✗



Results

ANOVA 2-way interactions


Source	Prob>F	Significant ($p \leq 0.05$)
Nodule Type	0	✓
NType X SMethod	0	✓
Sizing Method	0	✓
Readers	0	✓
Readers X SMethod	0	✓
Readers X NType	0.0002	✓



Results

ANOVA 2-way interactions

Source	Prob>F	Significant (p≤0.05)
Nodule Type	0	✓
NType X SMethod	0	✓
Sizing Method	0	✓
Readers	0	✗
Readers X SMethod	0	✗
Readers X NType	0.0002	✗




Results

ANOVA

- Important subgroups for analysis
 - Sizing method by nodule type
 - Compare sizing methods for
 - 10 mm Sphere
 - 20 mm Sphere
 - 20 mm Ellipsoid
 - 10 mm Lobulated
 - 10 mm Spiculated


Results

Bias comparisons




Results

Variance comparisons



Secondary analysis

Variance comparisons



- 0.8 mm slice data only

Part 2 Summary

- **Bias analysis**
 - Radiologist 3D volume estimates were close to unbiased across all nodule types
 - Radiologist 1D & 2D sizes systematically underestimated true size for complex nodules
 - Different biases for elliptical, lobulated and spiculated nodules
 - Not shown
 - Nodule orientation critical factor for 1D & 2D but not 3D
- **Variance analysis**
 - Variability lower for 1D/2D sizing compared with 3D volume
 - Not shown
 - Scaling 1D & 2D sizes to 3D resulted in more comparable variability across methods
- Each sizing method had its own unique bias/variance tradeoff

Acknowledgments

- We acknowledge Corelab's strong support to this QIBA vCT Part 1A groundwork effort. CoreLab Partners LLC conducted the reader study component of this investigation and CoreLab radiologists participated as readers
- We acknowledge the substantial contributions of Lisa M. Kinnard (Medical Research Program, Department of Defense, Fort Detrick, MD) in the design and conduct of the QIBA 1A reader study
-
- We acknowledge the members of the QIBA Volumetric CT Technical Committee and especially the members of the QIBA vCT Part 1A subcommittee for making substantial contributions to Part 3.

Acknowledgments

- The phantom data collection was funded through a Critical Path grant from the U.S. Food and Drug Administration. The intramural research program of the National Institute of Biomedical Imaging and Bioengineering and the National Cancer Institute through IAG no. 224-07-6030 provided partial support for the phantom data collection. Phantom scans collected on the Philips IDT Mx8000 were supported in part by the Center for Interventional Oncology at the National Institutes of Health (NIH) and an Interagency Agreement between the NIH and the United States Food and Drug Administration (FDA).
- The mention of commercial entities, or commercial products, their sources, or their use in connection with material reported herein is not to be construed as either an actual or implied endorsement of such entities or products by the Department of Health and Human Services.