

Image Archives

Where are they? What's next?

Lung Cancer Workshop

May 2 - 3, 2011

NIST Biochange:

Charles Fenimore, John Lu, Jim Filliben

Publicly available archives



- NBIA collections
 - Used in biochange benchmarking challenges
 - Growing with additions from FDA & QIBA
- IELCAP / Cornell / Mount Sinai
 - Volcano '09
- Give-A-Scan
- Potential – NLST imagery, NELSON trial, Pharma contributions, QIN

NBIA collections

- Growing with QIBA
- Public access
 - CT Colonography ACRIN
 - CT - LIDC (public) + IDRI (limited access)
 - NCRI (National Cancer Research Institute) – MR (NHS pilot)
 - Phantom - University of Washington Public
 - CT Phantom, FDA Public
 - QIBA DCE-MRI - ~5 Supplying Sites and 1 Analysis site. Limited
 - RIDER (MSKCC, MDACC, UWashington, Duke, UMichigan) (322 patients)
 - CT - Roswell Strong (ROSWELL-STRONG)
 - MR - TCGA Henry Ford, USCF
 - VASARI TJU, Henry Ford
 - Virtual Colonoscopy WRAMC, San Diego

Use of RIDER in Biochange Challenge

Single Slices from Validation CTs




New phantom imagery

QIBA vCT 1C


- Develop and apply two arms of the Imaging Protocol
 - ACRIN Study (6678)
 - Quality-based protocol: Aim is constant image quality across scanners
- Image anthropomorphic phantom at five imaging sites
 - 5 sites with 5 scanners from 4 manufacturers
 - Phantom to have 12 nodules with spherical and spiculated shapes
- Clinicians measure synthetic nodules
 - Semi-automated 3-D volume & derived 1- & 2-D sizing
- Analysis
 - Variability between scanners, arms of protocol & sizing measures
 - Inter-reader variation

QIBA vCT 1C

- Determining the image quality/performance-Based arm of protocol, device-independence



Target: 6 lp/mm



Target; $17 \pm 1\text{HU}$

QIBA 1C Groundwork

- Scanners

- GE VCT 64, Philips Brilliance 16, Philips Brilliance 64, Siemens Sensation 64, Toshiba Aquilion

- Status

- Anthropomorphic phantom data and Gammex CT phantom data collect on all scanners – 210 series
- Reader study under development

QIBA 1C sample imagery

