

QIBA 1A: Evaluation of nodule sizing by radiologists through CT phantom imaging of synthetic nodules

Nicholas Petrick, Hyun J. Grace Kim, David Clunie,
Kristin Borradaile, Robert Ford, Rongping Zeng,
Marios A. Gavrielides, Michael F. McNitt-Gray,
Charles Fenimore, Z. Q. John Lu, Binsheng Zhao
Andrew J. Buckler

Author Affiliations

Nicholas Petrick^a, Hyun J. Grace Kim^b, David Clunie^c, Kristin Borradaile^c,
Robert Ford^d, Rongping Zeng^a, Marios A. Gavrielides^a,
Michael F. McNitt-Gray^b, Charles Fenimore^e, Z. Q. John Lu^e, Binsheng
Zhao^f Andrew J. Buckler^g

^aCenter for Devices and Radiological Health, U.S. Food and Drug Administration

^bDavid Geffen School of Medicine at UCLA, Los Angeles, CA

^cCoreLab Partners, Inc, Princeton, NJ

^dPrinceton Radiology, Princeton, NJ

^eNational Institute of Standards and Technology, Gaithersburg, MD

^fDepartment of Radiology, Columbia University Medical Center, New York, NY

^gBuckler Biomedical LLC, Wenham, MA, USA

Outline

- Objectives
- Study Design
 - Dataset
 - Reading protocol
 - Analysis method
- Results
- Summary
- Ongoing QIBA Projects

QIBA vCT Groundwork

- Quantitative Imaging Biomarker Alliance (QIBA)
 - To investigate the role of quantitative imaging methods as potential biomarkers in evaluating disease and responses to treatment
 - Current Technical Committees
 - Volumetric CT (vCT)
 - COPD CT
 - DCE-MRI
 - fMRI
 - FDG PET-CT

Objectives

- Overall Objected of QIBA vCT Effort
 - To develop groundwork data evaluating quantitative CT tumor sizing methods as imaging biomarkers
 - Clinical Trails
 - Clinical Practice
- Current study is the QIBA vCT Part 1A groundwork effort
 - Objective
 - To estimate bias/variance of radiologists estimating the size of synthetic nodules from CT scans of an anthropomorphic phantom

Dataset (Nodules)


- 10 synthetic nodules
 - 5 shapes X 2 densities
- Shape
(Volume equivalent to sphere of given diameter)
 - 10 mm Sphere
 - 20 mm Sphere
 - 20 mm Ellipsoid
 - 10 mm Lobulated
 - 10 mm Spiculated
- Density
 - -10 HU
 - +100HU



Spherical



Ellipsoid




Lobulated



Spiculated

Dataset (Thorax Phantom)



Anthropomorphic thorax phantom (Kyotokagaku Incorporated, Tokyo, Japan)

Dataset (Nodule Attachment)

- Nodules attached to synthetic vasculature



Dataset (CT Scanning)

- CT Scanner
 - Philips 16-slice MxIDT 8000 scanner
- Acquisition Parameters
 - Exposure (120 kVp): 100 mAs/slice
 - Pitch: 1.2
 - Recon kernel: B60 (Detailed)
 - Slice thickness (50% overlap):
 - 0.8 mm (16X0.75 collimation)
 - 5.0 mm (16X1.5 collimation)
 - Repeat Exposures: 2 scans for each nodule

Dataset

- 40 total dataset for segmentation
 - 10 nodules X 2 slice thickness X 2 repeat scans

Reading Protocol


- Readers
 - 6 radiologists familiar with evaluating lesion response in drug trials
- Sizing Methods
 - 1D technique (linear distance)
 - Largest in-slice diameter for the lesion
 - Based on RECIST criteria
 - 2D technique (area)
 - Largest in-slice diameter for the lesion
 - Largest perpendicular diameter within same slice
 - Based on WHO criteria
 - 3D technique (volume)
 - Semi-automated volumetric measurement tool

Reading Protocol

- Semi-automated 3D lesion sizing
 - Define seed strokes
 - Apply segmentation tool
 - Evaluate quality of segmentation
 - Refine segmentation
 - Adding/subtracting seed strokes only
 - Repeat until satisfactory segmentation achieved
- Software provided volume estimate

Reading Protocol

Session 1



Session 2

- Same process
- Randomized reading order

Readers applied all 3 sizing techniques within each session
120 measurements per reader per session

Reading Protocol

- All readings session took place at CRO central facility
 - Proprietary software application
 - Consumer color LCD monitors
 - DICOM Grayscale Standard Display Function calibration
 - Lung window/level (1200 HU/-600 HU)
- Each measurement technique was independently applied (readers did not see their previous measurements concurrently)

Analysis

- Compared bias/variability among methods

$$Size_{Rel} = \frac{Size_{Est} - Size_{True}}{Size_{True}} \bullet 100\%$$

$$Bias_{Rel} = Bias(Size_{Rel})$$

$$Std_{Rel} = \frac{Std_{Est}}{Size_{True}}$$

Truth

- 1D
 - Longest diameter (calipers)
- 2D
 - Longest diameter (calipers)
 - Longest perpendicular diameter (calipers)
- 3D (weight-density method)
 - Measured weight (precision scale)
 - Density (from manufacturer)

Regression Model


- Mixed-effect linear regression
 - Model relative bias
 - Fixed Effects
 - Nodule shape/size nodule density, slice thickness
 - Random Effect
 - Readers
- Comparison of biases
 - 1D vs. 2D
 - 1D vs. 3D
 - 2D vs. 3D

Example Case


- 10 mm spiculated
 - True longest diameter
 - 22 mm
 - True area
 - 343 mm²
 - True volume
 - 529 mm³



Example Case




1D Size_{Rel} = -9.1%




2D Size_{Rel} = -44.6%

Example Case




3D In-slice



3D Size_{Rel}= +36.1%

Combined Results


Size Method	Relative Bias	Relative Std. Deviation
1D	-14.6%	20.4%
2D	-18.8%	28.3%
3D	-1.3%	21.9%



Comparison of Biases

- 1D vs. 2D
 - -14.6% vs. -18.8%, $p < 0.001$
- 3D vs. 1D
 - -1.3% vs. -14.6%, $p < 0.001$
- 3D vs. 2D
 - -1.3% vs. -18.8%, $p < 0.001$

- By shape
 - Spherical
 - Non-spherical
- And slice thickness
 - 0.8 mm
 - 5 mm



Bias and Variance

Size Method	Spherical Nodules		Non-spherical Nodules	
	0.8 mm	5.0 mm	0.8 mm	5.0 mm
1D	2% (± 5)	0% (± 4)	-23% (± 20)	-27% (± 21)
2D	4% (± 10)	0% (± 11)	-33% (± 26)	-33% (± 29)
3D	1% (± 12)	5% (± 23)	0% (± 14)	-2% (± 30)

Summary

- 3D volume provides a low bias estimate of nodule volume
 - For spherical and non-spherical nodules
- 1D and 2D low bias only for spherical nodules
- Variance analysis similar relative precision between the 1D & 3D
 - 2D method larger relative standard deviation
- **3D volume at 0.8 mm slices, low bias and low variance**
 - Independent of shape

Future Work

- Statistical analysis stratified by nodule and CT acquisition characteristics
- Statistical comparison of precision (variability)
 - Intra- and inter-reader
- Analysis of reader segmentation data
 - STAPLE or similar type of analysis

Acknowledgments

- We acknowledge Corelab's strong support to this QIBA vCT Part 1A groundwork effort. CoreLab Partners LLC conducted the reader study component of this investigation and CoreLab radiologists participated as readers
- We acknowledge the substantial contributions of Lisa M. Kinnard (Medical Research Program, Department of Defense, Fort Detrick, MD) in the design and conduct of the reader study
-
- We acknowledge the members of the QIBA Volumetric CT Technical Committee and especially the members of the QIBA vCT Part 1A subcommittee for making substantial contributions to this work.

Acknowledgments

- The phantom data collection was funded through a Critical Path grant from the U.S. Food and Drug Administration. The intramural research program of the National Institute of Biomedical Imaging and Bioengineering and the National Cancer Institute through IAG no. 224-07-6030 provided partial support for the phantom data collection. Phantom scans collected on the Philips IDT Mx8000 were supported in part by the Center for Interventional Oncology at the National Institutes of Health (NIH) and an Interagency Agreement between the NIH and the United States Food and Drug Administration (FDA).
- The mention of commercial entities, or commercial products, their sources, or their use in connection with material reported herein is not to be construed as either an actual or implied endorsement of such entities or products by the Department of Health and Human Services.