

Pulmonary Nodules: Challenging Cases

Anthony P. Reeves

Vision and Image Analysis Group

School of Electrical and Computer Engineering

Cornell University

Lung Cancer Workshop VIII, May 2, 2011

Outline


Nodule Challenges

1. Large complex lesions
2. Different scan resolutions (slice thickness)
3. Vessel attachments
4. Scanner artifacts/motion
5. Part-solid nodules
6. Size is an issue
7. Cardiac motion/respiration motion
8. Non-solid nodules

Easy Cases


Example 1 **Good** result Axial visualizations at three times

Axial View



Cornell University
Vision and Image Analysis Group

VIA



Cornell University
Vision and Image Analysis Group

VIA

Web-system growth tables for Case 1

Group 1: 3D automated analysis

Nodule Information

Nodule	Date	Size (mm) from X-Y Extent	Size (mm) from Volume	Volume (mm ³)	Resolution
T1	Jun 13 2000	9.38	8.63	336.19	0.1875 x 0.1875 x 1.00
T2	Jul 19 2000	10.75	9.19	406.25	0.1875 x 0.1875 x 1.00
T3	Sep 26 2000	11.63	10.48	603.27	0.1875 x 0.1875 x 1.00

Growth Information

Nodules	Date T1	Date T2	Interval (Days)	Volume			Notes
				% Change	DT	GI	
T1-T2	Jun 13 2000	Jul 19 2000	36	20.84%	131.78	17.1	
T1-T3	Jun 13 2000	Sep 26 2000	105	79.44%	124.47	18.2	
T2-T3	Jul 19 2000	Sep 26 2000	69	48.50%	120.96	18.8	



Cornell University
Vision and Image Analysis Group

VIA

Nasty Case


Spiculations and many vessels

T1-T2: 203 days


-

From the Prevent Cancer Foundation:
Public Lung Database To Address Drug Response

T1: Center + ROI



xres: 0.693 zres: 2.5 mm



60 mm


Selected Slices




Cornell University
Vision and Image Analysis Group

VIA

T1: Manual Marking




97 mm



Cornell University
Vision and Image Analysis Group

VIA

T1: Projections



60 mm

Coronal

60 mm

Sagittal



Cornell University
Vision and Image Analysis Group

VIA


Large Complex Nodule




Cornell University
Vision and Image Analysis Group

VIA

Complex T1: Central slice + ROI



xres: 0.546 zres: 2.5 mm




50 mm



Cornell University
Vision and Image Analysis Group


VIA

T1: Projections



50 mm

Coronal



50 mm


Sagittal




Cornell University
Vision and Image Analysis Group

VIA

T2: Central slice + ROI




xres: 0.187 zres: 1.0




50 mm

Different scan resolutions for T1 and T2: 0.546×2.5 vs. 0.187×1.0
Voxels in T1 are 21 times larger than voxels in T2

T1: Volumetric Segmentation



50 mm



Volumetric analysis worked for this case. It considered the central region (red) as solid and the surrounding area (blue) as sub-solid

Large Complex Nodules

- Standard volumetric methods designed for small nodules fail
- Solution: use a density based method that evaluates the change in density over a fixed size region as a surrogate for volume change
- Outcome: all nodules have successful size change estimates

Case	Interval (days)	% Size Change	GI (%/month)	DT (days)
0414	44	9.0	6.2	350
1297	56	21.2	11.0	202
2461	29	12.7	13.4	168



Cornell University
Vision and Image Analysis Group

VIA

Nasty Case


Zero-Change Dataset

T1-T2: ~5 minutes


-

From the VOLCANO '09 Challenge


T1 and T2 1.25 mm slice thickness



T1 and T2 1.25 mm slice thickness



Zero-change automated segmentation



Cornell University
Vision and Image Analysis Group

VIA

Scanner artifacts?

Issue: Different number of slices on scanner reconstructions

- Z-compensation method. A volumetric evaluation method that is insensitive to linear transformations in the Z dimension.
- Measured volume change:

Standard volumetric method	30.3%
Z-compensation method	12.6%




Cornell University
Vision and Image Analysis Group

VIA


Size Range

Size Issue: two different nodules



4.86 mm (60 mm³)

47 mm




29.14 mm (12961 mm³)

49 mm




Cornell University
Vision and Image Analysis Group


VIA



Segmentation (large)




Selected slices




6 cardinal views

Segmentation (small)



Selected slices

50 mm




6 cardinal views



Cornell University
Vision and Image Analysis Group


VIA

Size Issue



4.86 mm (60 mm^3)

$20 \text{ mm}^3 = 30\% \text{ change}$



49 mm

29.14 mm (12961 mm^3)

$20 \text{ mm}^3 = 0.14\% \text{ change}$

1:216



Cornell University
Vision and Image Analysis Group

VIA

Addressing Challenging Nodules

1. Issue: Very Complex Large Nodules
Use Density Change measurement method
2. Issue: Scanner slice variation (due to reconstruction artifacts?)
Use a method with Z-compensation
3. Issue: Very small Nodules
Use a detailed algorithm that corrects for differences in segmentation
4. Issue: Different scan parameters
Careful image resampling can help mitigate this issue
5. Issue: Patient/scanner motion, inconsistent data
Take additional scans



Cornell University
Vision and Image Analysis Group

VIA